


Sample Type


Sample Type allows users to configure new Sample Types and edit the existing ones.

Accessing Sample Type

- 1. Login into the application.
- 2. Select Quality.
- 3. Click on Sample Type.



- 4. Clicking on **New** display the following dialog



- 5. The fields present in this dialog are listed below:

Field	Description	Example
Sample Type	Enter the Sample Type	Tea
Description	Enter the Description of the Sample Type	Green Tea
Control Point	Select the Control Point	
Sample Attributes		

- 6. In order to edit a **Sample Type**, click the desired **Sample Type** in the **Sample Type** screen, to highlight it in the summary list.
- 7. All the available sample attributes are listed out in the Attribute Name drop down. The required attributes can be selected and deselecting using the **Remove** button.
- 8. Double click the highlighted row or click **View button**. The **Sample Type Detail** dialog opens.
- 9. The desired fields can be edited.
- 10. Clicking **Save** saves the Sample Type.
- 11. Clicking **Delete** deletes the Sample Type.