Test

Overview

- Test will be defined with multiple properties associated with it
 The user can create Test containing multiple Properties associated with it
 "Test" configuration screen can be accessed by clicking on "Test" option under "Quality"
 Test" configuration screen is shown below with sample test

Active
\checkmark
\checkmark

- Click on "New" to create new test
- Click on "View" to view / edit the selected test

Test Configuration

 $^{\circ}\;$ Test configuration requires the following parameters to be configured

Field	Description	
Test Name	Name of the test	
Description	Description of the test	
Analysis Type	Type of Analysis that requires for the test like Physical, Chemical, etc	
Test Method	Method for the Test	
	Ex.: Food Chemicals Codex 3rd Ed. method III	
Industry Standards	Industry Standards for the test	
Sensorial Label	Sensorial Label associated to the test	
Replications	No. of times the test to be conducted	
Auto Capture	Whether the test values will be auto captured or manual	
Active	Set to activate the test	
Ignore Sub Sample Count	When this value is set, the system will not make the "Replications" count mandatory	
Auto Capture	Whether the test values will be auto captured or manual	
Properties Association	Associating the relevant properties to the list	

🛟 Test - Tea Blend	
New Save Search Delete Undo Close	
Test Details	Properties
Test Name: Tea Blend	Remove H Layout Filter Records (F3) Property Name
Description: Tea Blend	PP-Test-Integer
Analysis Type: Physical	
Test Method: Tea Blend	
Industry Standards: 1	
Sensorial Label: 1	
Replications: 1	
Auto Capture:	
Active:	
Ignore Sub Sample Count: 📃	
🕝 😧 🖓 Ready	4 Page 1 of 1 ▶ ▶

Note: Test can not be deleted if it is already associated to a Quality Template