

# How do I create Test?

- 1. Navigate to Quality --> Test
- 2. Test Manager screen is loaded

Test

New View Refresh Close

Layout Filter Records (F3) 2 records

Test Name	Description	Test Method	Analysis Type	Industry Standards	Sensorial Label	Replications	Active
Tea Blend	Tea Blend	Tea Blend	Physical	1	1	1	<input checked="" type="checkbox"/>
Tea Receipt	Tea Receipt	Tea Receipt	Physical	Tea Receipt	1	1	<input checked="" type="checkbox"/>

- 3. Click on New
- 4. Test detail screen is loaded.

Test -

New Save Search Delete Undo Close

Test Details

Test Name:

Description:

Analysis Type: Physical

Test Method:

Industry Standards:

Sensorial Label:

Replications: 1

Auto Capture: ☐

Active: ☒

Ignore Sub Sample Count: ☐

Properties

Remove Layout Filter Records (F3)

Property Name

Edited Page 1 of 1

- 5. Enter the Test Name.
- 6. Enter the Description
- 7. Select the Analysis Type from drop down.
- 8. Enter the Test Method.
- 9. Enter Industry Standards
- 10. Enter the Sensorial Label.
- 11. Select the Properties from the Property Name drop down

Properties

Remove Layout Filter Records (F3)

Property Name

Name	Description	DataType Name	List Name
PP-Test-Integer	PP-Test	Integer	
PP-Test-Float	PP-Test-Float	Float	
PP-Test-Bit	PP-Test-Bit	Bit	
PP-Test-List	PP-Test-List	List	MouthFeel
PP-Test-String	PP-Test-String	String	
PP-Test-Time	PP-Test-Time	Time	
PP-Test-DateTime	PP-Test-DateTime	DateTime	

- 12. Click on Save.