


Contamination Group

This screen allows user to define contamination groups.

Contamination Group Creation:

- 1. Open Contamination Group
- 2. Clicking on **NEW** displays the following screen:



- 3. The fields in the Contamination Group screen are listed below:

Field Name	Description
Contamination Group	Enter the name of Contamination Group
Short Name	Enter short name for the group
Description	Enter description for the group

- 4. Clicking on **Insert** adds a new row in the grid.
- 5. Clicking on **Remove** removes the selected rows.
- 6. Clicking **Save** saves the information entered in the grid
- 7. Clicking **Undo** reverts the grid to last saved data.